 **ceramill motion****Bedienungsanleitung****User Manual****Mode d'emploi****Instruzioni d'uso****Modo de empleo**

Deutsch	2 - 23
English	24 - 45
Français	46 - 67
Italiano	68 - 89
Español	90 - 111

**AMANN GIRRBACH**

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Tab. 6

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<b>1</b>	<b>Explanation of Symbols</b> .....	<b>25</b>	<b>8</b>	<b>Cleaning and Maintenance</b> .....	<b>40</b>
<b>2</b>	<b>General Safety Instructions</b> .....	<b>26</b>	8.1	Cleaning .....	40
<b>3</b>	<b>Suitable Personnel</b> .....	<b>26</b>	8.2	Weekly Maintenance .....	40
<b>4</b>	<b>Machine Specifications</b> .....	<b>27</b>	8.3	Monthly Maintenance .....	43
4.1	Delivery Scope .....	27	8.4	External Maintenance .....	43
4.2	Intended Use .....	27			
4.3	CE Declaration of Conformity .....	27	<b>9</b>	<b>Malfunctions, Repairs and Warranty</b> ..	<b>44</b>
4.4	Components and Interfaces .....	28	9.1	Malfunctions .....	44
<b>5</b>	<b>Installation</b> .....	<b>30</b>	9.2	Repairs .....	44
5.1	Setting Up .....	30	9.3	Warranty .....	44
5.2	Service Unit .....	30			
5.3	Installation of the Ceramill Match and Ceramill Motion Software .....	32	<b>10</b>	<b>Environmental Protection</b> .....	<b>44</b>
5.4	Establishing the Network Connection .....	33	<b>11</b>	<b>Technical Data and Spare Parts</b> .....	<b>45</b>
5.5	Starting the Ceramill Motion Operation .....	33			
<b>6</b>	<b>Application and Operation</b> .....	<b>35</b>			
6.1	Checking the Model Area .....	35			
6.2	Function of the buttons .....	36			
6.3	Reference Run .....	36			
6.4	Program Execution .....	36			
6.5	Checking the Workpiece During a Running Program .....	37			
6.6	Moving the the Service Position .....	37			
6.7	Moving the the Tool-changing Position .....	37			
<b>7</b>	<b>Registration and Download Information</b> .....	<b>38</b>			
7.1	Quick Registration, M-center .....	38			
7.2	Order Tracking, M-Center .....	39			
7.3	Download Infos .....	39			

# 1 Explanation of Symbols

## Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- \_ **NOTE** means that property damage can occur.
- \_ **CAUTION** means that light to fairly serious personal injury can occur.
- \_ **WARNING** means that serious personal injury can occur.
- \_ **DANGER** means that serious personal injury can occur.

## Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

## Other symbols in the Manual

Symbol	Meaning
>	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Tab. 1

## Other symbols on the machine

Symbol	Meaning
	USB port (socket type B)
	Network connection (Ethernet)
	Slot for SD card
	Control cable connection, extraction
	Fuse
	Compressed-air connection

Tab. 2

## 2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

### CAUTION:

A defective machine can lead to malfunctions!  
On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

### NOTE:

Possible machine damage through insufficient extraction!

- ▷ Operate the machine only with Ceramill Air-stream or another compatible extraction system (see page 9).

### NOTE:

- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

## 3 Suitable Personnel

### NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

## 4 Machine Specifications

### 4.1 Delivery Scope

- \_ Ceramill Motion milling machine
- \_ Mains cable
- \_ USB cable (5 m)
- \_ Network cable (5 m)
- \_ Control cable, extraction
- \_ Compressed-air hoses (250 mm and 2000 mm)
- \_ Service unit
- \_ Service set for Jäger spindle
- \_ Ceramill-Match-Package
  - Ceramill Match CD (Ceramill Match and Ceramill Motion operating software, installation instructions)
  - Ceramill-Match-Dongle
  - User Manual
- \_ Blank template, ZI71 Motion
- \_ Ceramill TEST 71 L (test blank)
- \_ Ceramill Motion Roto 2.5 (milling tool)
- \_ Screw set
- \_ Blank holder
- \_ Extractor adapter

The Ceramill Motion is delivered with a Ceramill TEST 71L test blank (Order No. 760301) inserted in the blank holder. Prior to delivery, a test object was factory-milled in order to check the zero point settings. Possibly remaining milled shavings/chips are not a reason for rejection. The test blank is provided for test milling and calibration of the Ceramill Motion.

- ▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier (compare with Ceramill Mall setup checklist).

### 4.2 Intended Use

The Ceramill Motion is a PC-controlled milling machine for producing dentures by milling blanks of pre-sintered zirconium oxide, plastics and wax.

Using blanks and milling tools not approved by AmannGirrbach can damage the machine and make the milling product unusable. For such cases, AmannGirrbach shall assume no liability whatsoever.

Unauthorised modifications/alterations of the machine shall also void the warranty.

### 4.3 CE Declaration of Conformity

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

The product's Declaration of Conformity can be viewed under [www.amanngirrbach.com](http://www.amanngirrbach.com).

## 4.4 Components and Interfaces

EN


Fig. 1 Machine overview, front

- 1 Measuring device, tool length
- 2 Tool holder
- 3 Milling spindle
- 4 Blank holder
- 5 Dust tray
- 6 Service button
- 7 Operation button


Fig. 2 Machine overview, rear

- 1 USB port
- 2 Network connection
- 3 Programming interface RS 232
- 4 SD card slot
- 5 Socket for control cable, extraction (for automatic operation)
- 6 Compressed-air hose, 250 mm
- 7 Compressed-air inlet
- 8 Service unit
- 9 Compressed-air connection, milling spindle
- 10 Extraction-hose connection
- 11 Fuse compartment
- 12 Power supply connection
- 13 Main switch

## 5 Installation

### 5.1 Setting Up

- \_ The machine is intended exclusively for use within dry, closed rooms.
  - \_ Minimum space requirements for the Ceramill Motion (incl. connections):
    - Width: 615 mm
    - Depth: 430 mm
    - Height: 730 mm
  - \_ The set-up surface must withstand a load of 62 kg.
  - \_ During operation, the room temperature must be between 18 °C and 30 °C; avoid large variations in temperature.
- ▷ Set up the machine on a sturdy work table or work bench (off the ground and clear of walls).
  - ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.

#### Countries with 100-115 V mains voltage

When the mains voltage at the set-up location is 100-115 V:

- ▷ Pull out the fuse compartment ([11] in Fig. 2 on page 7) below the mains plug.
- ▷ Replace both T3.15 A fuses with the provided T6.3 A fuses.

### 5.2 Service Unit

The Jäger spindle is equipped with a sealing-air feature. This sealing air prevents shavings/chips and dust from entering the spindle. The service unit filters possible contamination in the sealing air and thus prevents costly machine damage.

## Installing the service unit


Fig. 3 Service unit

- 1 Compressed-air inlet
  - 2 Compressed-air hose, 250 mm
  - 3 Water separator
  - 4 Compressed-air connection, milling spindle
  - 5 Extraction-hose connection
- ▷ Mount the service unit on the rear side of the machine using the knurled screws.
  - ▷ Connect the 250 mm compressed-air hose to the right-hand connection (viewed from the front) of the service unit and to the compressed-air connection [4] of the milling machine.
  - ▷ Connect the 2000 mm compressed-air hose to the compressed-air inlet [1] of the service unit and the other end to the compressed-air supply of the laboratory.

## Compressor requirements

The compressor for the Ceramill Motion must meet the following minimum requirements:

- \_ Inlet flow rate: 100 l/min (3.54 CFM)
- \_ Compressed-air reservoir: min. 50 l (13.2 gallons)
- \_ Output pressure: 8 bar / 116 psi
- \_ Sound pressure level: max. 60 dB(A)

A suction capacity of 70 l/min (2.47 CFM) at 5 bar (72.5 psi) is required during operation.

## Requirements on air purity for the Jäger spindle

When another service unit than the one supplied is to be used, the following connection conditions according to DIN-ISO 8573-1 must be ensured:

- \_ Solid contaminants: Class 3 (particle size max. 5  $\mu\text{m}$ , particle density max. 5  $\text{mg}/\text{m}^3$ )
- \_ Water content: Class 4 (max. pressure dew point +3 °C, water content max. 6000  $\text{mg}/\text{m}^3$ )
- \_ Total oil content: Class 2 (max. oil content 0.1  $\text{mg}/\text{m}^3$ )

## Extraction requirements

When another extraction system than the Ceramill Airstream is to be used, it must meet the following requirements:

- \_ Suction capacity: 56.6 l/min
- \_ Suitability for zirconium-oxide dusts
- \_ HEPA microfilter (97.97 %), filter class H12, dust category M

## 5.3 Installation of the Ceramill Match and Ceramill Motion Software

**NOTE:**

Malfunction!

▷ Operate the milling machine only via PCs delivered by AmannGirrbach!

- ▷ Insert the Ceramill Match CD-ROM into the drive of the PC.  
The CD contains a PDF file with installation guide.
- ▷ Open the PDF document and follow the instructions in the guide.

### Connecting the interfaces


Fig. 4 Interfaces on the rear side of the milling machine

- 1 USB data port
  - 2 Ethernet interface
  - 3 Slot for SD memory card
  - 4 Control interface for extraction
  - 5 Main switch
  - 6 Power supply connection
  - 7 Extraction connection
- ▷ Make sure that the main switch [5] of the milling machine is in the 0 position.

- ▷ Check if the SD memory card is properly seated in the slot [3] on the rear side of the milling machine.
- ▷ Connect the control cable to the control interface for extraction [4] on the rear side of the milling machine and to the interface on the rear side of the extraction system (e.g. Ceramill Airstream).
- ▷ Insert the extraction hose of the extraction system being used (e.g. Ceramill Airstream) to the extraction connection [7] on the rear side of the milling machine.
- ▷ Connect the mains cable to power supply connection [6] of the milling machine as well as to a mains socket outlet.

**NOTE:**

Malfunctions caused by driver conflicts on the PC.

▷ Do not connect the Ceramill Motion milling machine and the Ceramill-Match-Dongle to the PC until after the software is installed.

- ▷ Connect a USB cable to the USB data port [1] of the milling machine and to a free USB port of the PC.  
The maximum permitted length of the USB cable is 5 m.

 The milling machine is only recognized by the software when it is connected to the PC via the same USB port, on which it was installed.

- ▷ Mark the USB port used for the milling machine on the PC, and use only this port for the connection.

- ▷ Insert the Ceramill-Match-Dongle in a free USB port on the rear side of the PC.


Fig. 5 Ceramill-Match-Dongle

### 5.4 Establishing the Network Connection

To establish the network connection, the milling machine and the PC must first be connected via USB.

- ▷ Connect the Ethernet interface [2] of the milling machine with the PC or with the network using a network cable.  
The maximum permitted length of the network cable to the PC or to the network router is 100 m.
- ▷ Switch the milling machine on.
- ▷ Start the Ceramill Motion Software.
- ▷ After initializing the milling machine in the software under *Settings > Communication*, set the connection type from USB to Ethernet and adapt the IP address, if required (see description of the Ceramill Motion Software).
- ▷ Switch the milling machine off and then on again.
- ▷ Initialize the milling machine in the software.
- ▷ After the software has detected the Ethernet connection, remove the USB cable.

When operating via a network, up to eight milling machines can be connected and controlled via the Ceramill Match Software.

### 5.5 Starting the Ceramill Motion Operation

In addition to standard tools, sister tools can optionally be inserted in the tool holder. When a standard tool has reached its wear limit, the Ceramill Motion automatically accesses the respective sister tool.

- ▷ Insert the Ceramill Motion Roto standard tools into the tool holder of the Ceramill Motion. Position the tool shank facing toward the outside.
  - Roto 2.5 is in the spindle.
  - Set Roto 1.0 to position T2.
- ▷ Optional: Insert the sister tools into the tool holder of the Ceramill Motion. Position the tool shank facing toward the outside.
  - Set Roto 2.5 to position T4.
  - Set Roto 1.0 to position T5.


Fig. 6 Tool holder

- ▷ Insert the selected blank into the holder in the Ceramill Motion. The high/thick side of the blank must face the spindle.


Fig. 7 Inserting a blank

- ▷ Fasten the blank using four screws and washers.

**NOTE:**

Do not damage the blank!

- ▷ Do not overtighten the screws! Avoid mechanical stress!

- ▷ Tighten the screws with the Allen key provided.


Fig. 8 Tightening a blank

- ▷ Shut the milling machine door.
- ▷ Switch the milling machine on via the main switch [4].
- ▷ Switch the PC on.
- ▷ Start the Ceramill Motion Software.  
When the software is started, a machine initialisation takes place.

If the milling machine is not detected:

- ▷ Restart the PC.

## 6 Application and Operation

### WARNING:

Danger from shavings/chips or tool fragments being thrown from the machine!

- ▷ Make sure that the door of the milling machine is always closed during the milling procedure!

### NOTE:

Insufficient milling results and increased wear!

- ▷ Do not operate the milling machine without extraction (e.g. Ceramill Airstream)!

The instructions for use (video tutorial) of the Ceramill Match Software and the User Manual for operation of the machine are on the Ceramill-Match-CD or can be downloaded under [www.ceramill-m-center.com](http://www.ceramill-m-center.com).

### 6.1 Checking the Model Area

The model area, on which the work is to be machined, must fit within the contour of the enclosed template (also see match package). Otherwise, the work cannot be machined in the blank. The template takes into account the shrinkage of the zirconium oxide in the final sintering.


Fig. 9 Template for ZrO<sub>2</sub> zirconium oxide

- ▷ For work in wax or plastic, measure out the model area.

**6.2 Function of the buttons**


Fig. 10 Buttons

Button	Door shut	Door open
Operation button [A]	Starts or stops the reference run; starts or stops the current program; activates the control position	Activates the extraction
Service button [B]	activates the service position; activates the tool-changing position	Activates the sealing air

Tab. 3

 Instead of with the buttons, these functions can also be actuated via the software (see instructions of the Ceramill Motion Software).

**6.3 Reference Run**

After switching on, a reference run must be carried out to determine the mechanical zero points. Operation button [A] flashes.

To start the reference run:

- ▷ Briefly push operation button [A].

The reference run is started. Operation button [A] flashes.

The reference run can be stopped:

- ▷ Briefly push operation button [A].

The reference run is stopped. Operation button [A] flashes.

To continue the reference run:

- ▷ Briefly push operation button [A].

The reference run is continued at the position where it was stopped. Operation button [A] flashes.

On completion of the reference run, the machine returns to the starting position. Operation button [A] goes out.

**6.4 Program Execution**

Before executing a program, it must first be loaded (see Ceramill Match Software). A loaded program is indicated by the lit operation button.

- ▷ Briefly push operation button [A].

The program is started. Operation button [A] lights up.

The program can be stopped:

- ▷ Briefly push operation button [A].

**-or-**

- ▷ Open the door.

The program is stopped. Operation button [A] flashes.

To continue the program:

- ▷ Briefly push operation button [A].

**-or-**

- ▷ Shut the door and briefly push operation button [A].

The program is continued at the position where it was stopped. Operation button [A] lights up.

On completion of the program, the machine returns to the starting position. Operation button [A] goes out.

## 6.5 Checking the Workpiece During a Running Program

▷ Push operation button [A] for at least three seconds.

The program is stopped. The spindle returns to the rear and the workpiece is turned around.

While travelling, operation button [A] is lit; after reaching the position, it flashes.

▷ Visually check the workpiece.

▷ Briefly push operation button [A].

The spindle travels to the position, where the program was interrupted. While travelling, operation button [A] is lit; after reaching the position, it flashes.

▷ Briefly push operation button [A].

The program is continued. Operation button [A] lights up.

## 6.6 Moving the the Service Position

The service position can be moved to from the starting position or from a stopped program.

▷ Briefly push service button [B].

The current tool is placed down and the spindle returns to the service position. The spindle is released; the motor is blocked. While travelling, operation button [A] is lit; after reaching the position, it flashes.

To exit the service position:

▷ Briefly push operation button [A].

The tool is picked up again. The spindle travels back to the starting position or to the position, where the program was interrupted. While travelling, operation button [A] is lit; after reaching the position, it flashes.

## 6.7 Moving the the Tool-changing Position

The tool-changing position can be moved to from the starting position or from a stopped program.

▷ Push service button [A] for at least three seconds.

The current tool is placed down and the tool holder is turned. While travelling, operation button [A] is lit; after reaching the position, it flashes.

To exit the tool-changing position:




▷ Briefly push operation button [A].

The tool is picked up again. The spindle travels back to the starting position or to the position, where the program was interrupted. While travelling, operation button [A] is lit; after reaching the position, it flashes.

## 7 Registration and Download Information


### 7.1 Quick Registration, M-center

Before registering, it is advisable to note down the serial numbers of the various equipment/devices.

	Device	No. location	Example	Serial number
[1]	Scanner (Ceramill Map)	Rear side of the unit		
[2]	PC	Top of PC housing		
[3]	Mind Dongle	Software box		
[4]	Motion (Ceramill Motion)	Rear side of the unit		
[5]	Match Dongle	Software box		

Tab. 4

- ▷ To register, log on to [www.ceramill-m-center.com](http://www.ceramill-m-center.com). On this page, click on the “registration” link. Afterwards, fill out the form field with the personal data.
- ▷ Enter all serial numbers.

Manufacturer: AmannGirrbach

Serial number of scanner:

Serial number of PC:

Serial number of Mind Dongle:

Serial number of Motion:

Serial number of Match Dongle:

- ▷ Accept the terms of use.
- ▷ Click on the “Create account now” button. A confirmation E-mail will be generated after a short period and sent to the previously entered E-mail address.
- ▷ Open the E-mail and click on the link in it. This opens a Web page that confirms the registration.

 If the link cannot be clicked on:

- ▷ Copy and paste the link into your Web browser.

Your e-mail address and a password of your choice are the access data for entering the protected homepage site [www.ceramill-m-center.com](http://www.ceramill-m-center.com) (order tracking, video tutorials).

Within the next 24 hours you will receive another E-mail. It contains the registration confirmation and the FTP access data (6-digit) for uploading data from our M-Center.

Two sets of access data are provided after the registration. One for the M-Center website (first E-mail) and one for the upload of your data to our M-Center (second E-mail).

▷ Enter all access data into the Ceramill Mind design software (the window opens automatically after clicking on the Upload button ).

### 7.3 Download Infos

Further instructions (video tutorials) on operating the Ceramill Map scanner and the Ceramill Mind software are available on the Instructions DVD or can be downloaded under [www.ceramill-m-center.com](http://www.ceramill-m-center.com).

Software updates are also available on our website under [www.ceramill-m-center.com](http://www.ceramill-m-center.com).

The registration process is now finished.

### 7.2 Order Tracking, M-Center

After receiving the design data in the M-Center, an order confirmation will be sent via E-mail.

Additionally, you will be informed via E-mail when your work goes into production and is dispatched.

A detailed list of all orders will be displayed on the M-Center homepage in the login area under "My orders". Here, each order will be allocated a tracking number after the order has been dispatched, with which it can be monitored.

Order number:	Completed
Name:	Ceramill 25L
Received:	10.05.2011 17:26:55
Planned date of dispatch:	12.05.2011
Cost:	€12.100%
STL file:	2010-11-22_00512-000-07-m-center_jed.stl
Stk. Trackingnumber:	Stamark number: 00112231010

## 8 Cleaning and Maintenance

### 8.1 Cleaning

The milling machine must be cleaned after each working day.

In order to clean the milling chamber and the spindle:

- ▷ Open the milling machine door.
- ▷ Push operation button [A] on the milling machine.

The extraction system is activated and the operation button lights up.

- ▷ Push service button [B] on the milling machine.

The sealing air at the spindle is activated and the service button lights up.

- ▷ Remove the dust tray and dispose of accumulated milling dust/shavings away from the machine.

In this manner, the filter bag in the extraction system will not fill up as quickly and can be used longer.

- ▷ Remove the extraction hose at the rear side of the machine and vacuum clean the milling chamber.

### 8.2 Weekly Maintenance

**NOTE:**

Danger of damaging the machine!

- ▷ Do not use compressed air, ultrasonic or a steam jet for the cleaning!

#### Spindle

The provided Jäger spindle service set and a milling tool of your choice are required for maintenance of the spindle.


Fig. 11 Service set for Jäger spindle

- 1 Collet grease
- 2 Collet brush
- 3 Felt cone
- 4 Ejector pin
- 5 Spindle wrench

In order to clean the spindle:

- ▷ Briefly push service button [B].  
The current tool is placed down and the spindle returns to the service position. The spindle is released; the motor is blocked.
- ▷ Open the milling machine door.
- ▷ Push service button [B] on the milling machine.  
The sealing air is activated and the service button lights up.
- ▷ Loosen the Allen screw at the bottom of the spindle cap and remove the cap.


Fig. 12

▷


Fig. 13

- ▷ Remove any dust from the spindle using a dry, soft cloth. Do not use cleaning agents!

To remove the collet:

- ▷ Apply the black spindle wrench [5] over the collet.


Fig. 14

- ▷ Insert a milling tool into the open collet.
- ▷ Using the spindle wrench, unscrew the collet out of the spindle.


Fig. 15

- ▷ Clean the collet holder using the felt cone [3].  
The collet holder must be free of shavings/chips and debris.


Fig. 16

- ▷ Clean the inside of the collet with the fine collet brush [2].


Fig. 17

- ▷ Apply a light coat of the collet grease ([1] in Fig. 11 on page 18) from the service set to the outer cone of the collet.  
This improves the sliding properties and increases the clamping force of the collet.

**NOTE:**  
Only use the provided collet grease for lubrication of the collet cone. Do not use other types of grease.

- ▷ Screw the cleaned and lubricated collet in reverse order to the stop back into the spindle.
- ▷ Mount the spindle cap again with the screw at the bottom.

**CAUTION:**  
Danger of tool breakage!  
▷ Make sure to remove the milling tool from the collet chuck!

When the milling tool was from the Ceramill Motion tool holder:

- ▷ Place the tool properly positioned back into the tool holder.
- ▷ Shut the milling machine door.
- ▷ Briefly push operation button [A].  
The tool is picked up again. The spindle travels back to the starting position or to the position, where the program was interrupted.

### Service unit

- ▷ Visually check the service unit for:
  - Visible particles, debris or particulate matter within the acrylic glass (signs for contaminated compressed air)
  - Amber discolouration of the filter element (sign of oil in the compressed air)
  - Water in the filter (sign of water in the compressed air)

Once any of the listed impairments is given:

- ▷ Correct the problem with the compressed air.
- ▷ Replace the complete service unit.

### 8.3 Monthly Maintenance

#### Checking the zero point

It is recommended to check the zero point once every month.

▷ Mill a test cube from the Ceramill TEST 71L test blank.

The test blank is to be used exclusively for test milling and calibration of the Ceramill Motion. The corresponding instructions are on the Ceramill Match CD or available on our website under [www.ceramill-m-center.com](http://www.ceramill-m-center.com).

#### Jäger spindle

The spindle must be in the starting position.

▷ Manually turn the shaft of the spindle through (at least ten times).


Fig. 18 Turning the Jäger spindle through

For further information on the Jäger spindle, please download the Manual from our website under [www.amanngirrbach.com](http://www.amanngirrbach.com).

### 8.4 External Maintenance

External maintenance is necessary after 1500 operating hours. The Ceramill Motion Software will inform you appropriately.

▷ Contacting the Ceramill helpdesk (workdays from 8 to 5 p.m.).

- Germany: +49 7231 957 100
- Austria: +43 5523 62333 390
- International: +43 5523 62333 399
- [helpdesk@amanngirrbach.com](mailto:helpdesk@amanngirrbach.com)

## 9 Malfunctions, Repairs and Warranty

### 9.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the milling machine.
- ▷ Restart the PC.

If these measures do not correct the problem:

- ▷ Contacting the Ceramill helpdesk (workdays from 8 to 5 p.m.).
  - Germany: +49 7231 957 100
  - Austria: +43 5523 62333 390
  - International: +43 5523 62333 399
  - helpdesk@amanngirrbach.com

### 9.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

### 9.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

## 10 Environmental Protection

### Packaging

In terms of packaging, AmannGirrbach participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

### Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the milling machine via public the public disposal systems in an environmentally-friendly manner.

The plastic parts are labelled. Thus, they can be sorted for recycling or disposal.

## 11 Technical Data and Spare Parts

 Subject to changes.

### Technical Data - Ceramill Motion

	Unit	Ceramill Motion
Art. No.		179200
Dimensions (D x W x H)	mm	465 x 430 x 730
Weight	kg	62
Electrical connections	V/A	230/3.15 100-115/6.3
Power output	W	250
Accuracy	µm	< 10
Allowable temperature range (operation)	°C	18 - 30
Spindle torque	Ncm	4
Spindle speed	rpm	60,000
Collet	mm	3
Sound pressure level	dB(A)	60
Compressed-air connection	-	Dry, clean compressed air
_ Pressure	bar	7
_ Flow rate	l/min	max. 70
Axes		4
Interfaces		USB / Ethernet / slot for SD card / RS 232 / Control cable for extraction system

Tab. 5

### Spare parts

Art. No.:	Designation
179210	Service unit, Motion
179211	Screw set for blank holder, Motion
179216	Threaded inserts for blank holder, Motion
179217	Test specimen for Motion test cube
179218	Service set, Jäger spindle
179229	Tool holder cover, Motion
179230	Collet 3 mm, Motion

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**Amann Gırrbach AG**
**Amann Gırrbach GmbH**