

# Powder X-Ray Diffraction Alignment Procedure

---

Procedure for aligning Rigaku Ultima III w/ Mercury CCD using  
Scintillation Detector for calibration

© 2013-14 **Thayumanasamy Somasundaram**

**11/18/2014**

This document describes the alignment procedure for Florida State University Chemical Sciences Lab's Rigaku Ultima III powder x-ray diffractometer with Mercury CCD. The procedure involves swapping the CCD with Scintillation Counter for alignment and back again for nano material powder x-ray diffraction. © 2013-14 Thayumanasamy Somasundaram.

## CONTENTS

Ultima III Alignment Procedure v0.9.6 (Beta) .....	2
Switching from CCD to scintillation counter and back .....	2
Start the shutdown .....	2
Swap the equipment .....	5
Swap CCD w/ Scintillation Counter and Micro Area w/ Standard Sample Holder .....	5
Restart .....	11
Restart w/ Scintillation Detector for alignment .....	11
Start the alignment procedure .....	13
Alignment procedure using .htm file .....	13
Complete Alignment, Save & Print Results .....	22
Swap the equipment back to original condition .....	24
Summary .....	27

© 2013-14 Thayumanasamy Somasundaram | Kasha Laboratory

Version 0.9.6 dated November 18, 2014

Institute of Molecular Biophysics | Florida State University, Tallahassee, FL 32306-4380

[tsomasundaram@fsu.edu](mailto:tsomasundaram@fsu.edu) | [www.biophysics.fsu.edu/soma](http://www.biophysics.fsu.edu/soma) | 1-850-644-6448

## ULTIMA III ALIGNMENT PROCEDURE V0.9.6 (Beta)

VER 0.9.5; DATED NOVEMBER 18, 2014

### SWITCHING FROM CCD TO SCINTILLATION COUNTER AND BACK

Written by Dr [Thayumanasamy Somasundaram](#) on November 17, 2004

The following write-up with photographs will show the alignment procedure to be done on Rigaku Ultima III Powder X-Ray Diffractometer (ru-pxrd) once in couple of months. Thanks to David Carnevale for help in demonstrating the procedure.

**Note:** This write-up is a beta version and some refinement is needed on some steps. Steps done at the PC (computer) will be shown in Times New Roman font; Sections done at the instrument r-pxrd will be shown in Georgia font.

### START THE SHUTDOWN

Rigaku Ultima III Powder X-Ray Diffractometer (shown on the right hand panel) with Mercury CCD and Scintillation Detector is located in CSL 1011 in FSU's Chemical Sciences Laboratory.

The generator has a Copper fixed anode tube and is cooled by Haskris chiller (originally a refrigerated version; now water-to-water chiller). The Mercury CCD is cooled by Julabo MD F25 chiller.



Haskris chiller for X-Ray tube



Julabo MD F25 chiller for CCD



Ultima III Powder XRD  
w/ Mercury CCD

Two computers control data collection and processing: 1) Control PC on the right and 2) Frame Grabber PC on the left. However, there is only ONE monitor and key board. So you access each computer using a KVM switch.

In the keyboard press, "Scroll Lock" key twice quickly to switch between the two PCs. Almost all of the time the monitor outputs and keyboard inputs are from and to Control PC.

Log in to control PC. It is likely that it is already running.



Monitor

Frame  
Grabber PC


Control PC

Keyboard

In the 'Standard Measurement Window', select

- File > Save
- Confirm

This allows the Standard Measurement values to be remembered when we restart the computer.




Close RaxVideo window



Left Click on RGKExpand on the left side of System Tray, then select


- Stop AreaMax

This closes the AreaMax window (better way to close it)




Close Standard Measurement Window

Close Manual Measurement Window



Left Click on RGKExpand on the left side of System Tray, then select

- Exit



In the Ultima III Generator's Left Top Panel, under the SAFETY RELEASE panel, press the 'DOOR' button. While the 'DOOR' button is blinking, open the sliding door all the way to left and right.




Inside the x-ray enclosure, locate the DALSA power supply (usually located on the right hand side and sits on the base plate) for the CCD detector. Power it down by switching it OFF (the power switch is located in the front of the unit, however, only the unit's back is visible) by locating the switch on the far right hand side.




Now, in the computer using the single keyboard press, "Scroll Lock" key twice quickly to switch to Frame Grabber PC.

- Log in using username and password
- Wait for the computer to display a DOS window
- Using the Start Button, shutdown the computer



Now, in the Ultima III Generator's Left Top Panel,


- Under X-RAYS
  - Press the Grey OFF button, then
- Under POWER
  - Press Red OFF button



Then, in the Ultima III Generator's Left Bottom Panel, open the door and locate the SERVO AMP UNIT,

- Turn it OFF by pressing the black power switch to off position.

Turn-off the Control PC now.



## SWAP THE EQUIPMENT

### SWAP CCD W/ SCINTILLATION COUNTER AND MICRO AREA W/ STANDARD SAMPLE HOLDER

With Ultima III generator door open, locate the CN ATT2 cable (this is plugged underneath the stage) and unplug it





Attach the CN ATT2 cable to the ATTACHMENT ADAPTER sitting on the base plate.

*The CN-ATT2 cable may fit loosely but make sure it secured as firmly as practical.*




Locate the plastic lid for the CCD and cover the CCD with it



Remove the micro-area stage (used for powder samples)

- Use the small hex wrench and undo the two screws until they are flush (on-top) with surface [2 or 3 mm hex wrench]
- Gently pull the stage toward the user
- Sit it down on the base plate of the instrument




Stage sitting on base plate

Remove the CCD Camera (with lid on it)


- Hold the CCD in one hand
- Note down the value on the back side (~105mm)
- Use Philips screw driver and undo all the four screws flush to the top
- Slowly wiggle the CCD and remove it
- Gently place it on the base plate



Stage and CCD sitting on base plate



Using larger hex wrench undo all the larger screws flush to the top and remove the black bracket that was holding the CCD and remove it completely (with the screws) and leave it on top of the CCD inside the enclosure.



Now, remove the two black brackets from the cabinet above the Control PC (they are marked #1 and #2)

- Put the #1 bracket close to the sample
- Put the #2 bracket away from the sample
- Use the hex wrench and secure them to the detector-arm of the instrument (the same place where the long black bracket was secured)



Now, we to put the Filter Holder (Part #53D-0466) in front and Scintillation Counter (Part # 57 38E401) and at the back with very little gap between them.



First, put the Filter Holder (FH) in front

- Lightly secure it in the detector-arm with 3 small hex screws. DO NOT fully tighten it
- Look underneath the FH and locate a "Red Line"
- Match this "Red Line" with the "Red Line" underneath the detector-arm
- Now tighten the Filter Holder fully



View from the bottom of the Filter Holder

Filter Holder in position now



Now, put the Scintillation Counter behind the Filter Holder with very little space between them

- Secure it with 2 hex screws

Now we have both Filter Holder and Scintillation counter in place



Plug in the connector for Scintillation counter



Plug in the RJ1 underneath the Filter Holder



Now, put the new stage Standard Sample Holder (stored inside the cabinet above the PC)


- Push the stage into the slot where the Micro-Area stage was installed
- The correct position of stage is guaranteed if the three holes in the stage are aligned with three holes on the instrument back panel
- Put the thumb screws and gently tighten
- Insert center slit

Center Slit and bubble level




With the center slit and bubble level in place now gently rotate the stage until the bubble is level both horizontally and vertically. Then hand tighten the thumb screws. The standard sample holder is now set.

- Remove the level and the center slit



Now open the Right Hand Side Bottom panel on the generator

- Change the SHT MODE switch from CCD to SD (flit it down)



Now open the Left Hand Side Bottom panel on the generator

- Press the first switch from bottom on the left-most PC card once.
- Briefly the blinking “green” LED light will stop blinking confirming your reset
- You will also see all the front panel LEDs light up and go away




- In CBO mirror box remove the H2 slit and put H10 slit. The lable should be facing the sample.
- Leave the PB in the holder



In the back filter holder

- Open the box
- Put Al foil in the front
- (0.3) slit in the back
- Put in 5 deg block (letter facing the user)
- Ni foil at the back
- All labels facing the sample

Close the box





# RESTART

## RESTART W/ SCINTILLATION DETECTOR FOR ALIGNMENT

Restart in this sequence


1. Switch ON the SERVO AMP UNIT
2. Press GREEN button in POWER Panel
3. Press White button on X-RAY panel
4. **Confirm Shutter 2 is in EXT mode**

Turn on the Control PC


Hit 'F1' to continue

Log in


- Administrator
- Password
- Right Click on Blue Rigaku Measurement server in Toolbar Menu
- Select RINT2200 and select to 'ShowWindow'
- Select Right and select to 'ShowWindow'



RINT2200 Window




Right Window




From Desktop open the Rigaku Control Panel Folder. In Control Panel open "System Details Button"

System Details Property Window change the XG Server Mode: **from** "Standard (Shutter Relay)




**To** "Standard"



Click 'Ok' and 'Confirm' the change. Right window will update.

Now click on the RINT2200 Right System (but not the ~~RINT2200 Right System Add/Delete button~~). This should open a new dialog box. In the following three tabs RINT2200 Right System Property need to be changed.

- System Construction
- X-Ray Beam Type
- Geometry System



**System Construction**

CHANGE FROM	CHANGE TO
<i>Goniometer:</i> Ultima 3 theta-theta goniometer+Mercury CCD	<i>Goniometer:</i> Ultima 3 theta-theta goniometer
<i>Attachement:</i> Micro area	<i>Attachment:</i> Standard Sample Holder
<i>Filter:</i> Not installed	<i>Filter:</i> K-beta filter
<i>Detector:</i> Mercury CCD	<i>Detector:</i> Scintillation Counter
<i>Temp Contoller:</i> Not installed	<i>Temp Contoller:</i> Not installed
<i>Incident Monochromateor:</i> use	<i>Incident Monochromateor:</i> use
<i>Slit:</i> Auto Variable (Mercury CCD/general)	<i>Slit:</i> Auto Variable Slit (Parallel/Ultima3)
<i>Monochromator:</i> Not installed	<i>Monochromator:</i> Not installed


**X-Ray Beam Type**

CHANGE FROM	CHANGE TO
No change. <i>Focus Selection:</i> Line; <i>Wave length selection:</i> K-alpha	



**Geometry System**

CHANGE FROM	CHANGE TO
<i>Geometry Information:</i> Micro Area	<i>Geometry Information:</i> Parallel Beam

After the change the three tabs should like the following







## START THE ALIGNMENT PROCEDURE

### ALIGNMENT PROCEDURE USING. HTM FILE

Now open the Auto-align(Right system).htm link in the Desktop


- It will open an Internet Explorer window and will load lots of stuff
- In the window, select the following
  - 2-Theta Alignment
  - Theta Alignment
  - Profile Measurement



Automatic alignment - Microsoft Internet Explorer

File Edit View Favorites Tools Help

Go Bookmarks 63 blocked Check AutoLink AutoFill Settings



Intensity (cps)

1000  
800  
600  
400  
200  
0

0.000 50.000 100.000 150.000

2 theta (deg.)

Grid Print Save Result Process Help Execute Counting Loss

2Theta Alignment  
 Theta Alignment  
 HV/PHA Adjustment  
 Monochromator Adjustment  
 Profile Measurement

Voltage 1: 40 kV Current 1: 40 mA HV/PHA Mono Sample d: 3.13532 A

Voltage 2: 40 kV Current 2: 40 mA Monochro Sample d: 3.3756 A

HV: 694.7 V Profile Sample d: 1.91999 A

XG ending  
 Present  
 Aging  
 X-ray off  
 Power off

Comment

Goniometer : Ultima3 theta-theta goniometer  
 Attachment : Standard sample holder  
 Slit : Auto Variable slit(Parallel/Ultima3)  
 Counter : Scintillation counter  
 Monochromator : Not installed  
 X-ray tube : Cu (Line focus)

Zero Position: -1.1750 mm

Zero Position :	-1.1750 mm	PHA Window :	80 div
Theta-D :	-0.2031 deg.	Monochro 2Theta :	--- deg.
Theta-S :	1.1104 deg.	Monochro Theta :	--- deg.
Theta-X :	--- mm	Monochro ThetaX :	--- mm
HV :	695 V	tIntegral :	0.7220 micro s
PHA Baseline :	60 div	tDifferential :	0.1833 micro s

Done My Computer

Now click the button “Execute” in the .htm file (center area) and follow the directions

- Confirm everything
- Initializing ... .. (takes sometime)
- X-Ray tube (left side of the instrument) will move
- 2 Theta arm (right side of the instrument) will move
- CBO box will move up and down
- 20 kv and 2mA will go to 40 kv and 40 mA
- Tube and 2 Theta will move to zero positions
- A scan of  $-1.0^{\circ}$  to  $+1.0^{\circ}$ Theta-D will be executed and the results will be displayed on the top portion of the htm file and bottom left hand side will also update information constantly

Automatic alignment - Microsoft Internet Explorer

File Edit View Favorites Tools Help

Google G Go 63 blocked Check AutoLink AutoFill Settings

Intensity (cps)

25000  
20000  
15000  
10000  
5000  
0

-1.000 -0.500 0.000 0.500 1.000

Theta-D (deg.)

Grid Print Save Result Process Help Execute Counting Loss

2Theta Alignment Voltage 1: 40 kV Current 1: 40 mA HV/PHA Mono Sample d: 3.13532 A  
 Theta Alignment Voltage 2: 40 kV Current 2: 40 mA Monochro Sample d: 3.3756 A  
 HW/PHA Adjustment HV: 694.7 V Comment Profile Sample d: 1.91999 A  
 Monochromator Adjustment  
 Profile Measurement

FXG ending  
 Present  
 Aging  
 Xray off  
 Power off

Goniometer : Ultima3 theta-theta goniometer  
Attachment : Standard sample holder  
Slit : Auto Variable slit(Parallel/Ultima3)  
Counter : Scintillation counter  
Monochromator : Not installed  
X-ray tube : Cu (Line focus)

Zero Position:-1.1750 mm

Zero Position :	-1.5000 mm	PHA Window :	100 div
Theta-D :	0.0000 deg.	Monochro 2Theta :	... deg.
Theta-S :	0.0000 deg.	Monochro Theta :	... deg.
Theta-X :	... mm	Monochro Theta'X :	... mm
HV :	695 V	Integral :	0.7220 micro s
PHA Baseline :	50 div	Differential :	0.1833 micro s

Done My Computer

**VBScript: Alignment confirmation**

Alignment a Multilayer mirror position.

Confirm the x-ray intensity indicated in the screen.

Operate mirror controller and find the position where the X-ray intensity hits maximum.

Click the <Stop> button when the x-ray intensity reaches the maximum.

Click the <OK> button to open the shutter.


It shows the intensity in the measuring progress in red in the profile frame.

OK Cancel

Now Press "OK" to confirm. The htm window will update and will now proceed to next step.


Now pick-up the DC Motor Controller (usually found on top of the instrument cabinet, on the Julabo side)

- Put two 9V batteries on the side opening
- Switch the Power ON (top of the controller)
- Now toggle the joystick left or/and right while watching either the .htm window or the right window
- The goal is to maximize the intensity (the joystick has some hysteresis and will occasionally reduce the intensity precipitously even for a small L or R toggle). So be patient and careful!
- Once successful, hit “Break” stop measurement and Confirm “Stop”
- The results are shown below.




Now follow the on screen instructions and do the following:

- Remove the Ni foil (K-beta filter); Open the door, remove filter, **close door of the filter box too!**
- Hit “Okay” on the pop-up window
- New htm window will open and three lines will be shown:
  - Red: Current Measurement (w/o K-beta filter!)
  - Green: Value with K-beta filter
  - Blue: Target intensity
- Once again the toggle the joystick to maximize the intensity. As close as possible to the “Blue line”.



Once maximum is achieved, click the “Break” button which stop the optimization procedure. Confirm to stop measurement. Typical result is shown below




Since no more optimization is possible, select “No More alignment” button when presented. Now switch “OFF” the DC Motor Controller, remove the two 9-VDC batteries and store the batteries in the cabinet box and place the controller back on top of the X-Ray Safety Cabinet (on the Julabo chiller side)





Center slit



Ni filter in the filter box

Now put the center slit in the standard sample holder and put the Ni-filter (aka K-beta filter) in the filter holder box. Click "OK" on the pop-up window to allow the measurement to continue. The first measurement result will be on Zero Position and is shown below. **There may be a Theta S Measurement after this!!!**

The screenshot shows the 'Automatic alignment' software interface. The main plot displays Intensity (cps) on the y-axis (0 to 8000) and Theta-S (deg) on the x-axis (-1.000 to 2.000). A red curve shows a sharp peak at 0.000 degrees. A 'Right Console' pop-up window is visible, displaying a radiation warning icon and the message: 'Theta-S Now scanning -1.0000 deg. to 2.9600 deg. !!'. Below the plot are various control buttons and settings.

**Right Console Message:**

```

Theta-S
Now scanning -1.0000 deg. to 2.9600 deg. !!
  
```

**Software Settings:**

- Voltage 1: 40 kV, Current 1: 40 mA, HV/PHA Mono Sample d: 3.13532 Å
- Voltage 2: 40 kV, Current 2: 40 mA, Monochro Sample d: 3.3756 Å
- HV: 694.7 V, Profile Sample d: 1.91999 Å


**Log Output:**

```

Execute setting HV. (695 V)
Execute setting PHA Base. (50 div)
Execute setting PHA Window. (100 div)
Execute setting XG power. (40 kV, 40 mA)
Execute moving Theta-S. (1.0000 deg)
Execute measuring Theta-D.
(START=-1,STOP=1,STEP=0.02,SPEED=4)
Execute moving Theta-D. (-0.1000 deg)
Rs. Slit = OPEN
Execute setting RSlit.
Measuring the X-ray intensity.
Measuring the X-ray intensity.
126582 cps
Measuring the X-ray intensity.
Execute measuring Ts.
(START=-2,STOP=2,STEP=0.01,SPEED=4)
Ts Correction:-0.350 mm
Execute moving Ts. (-0.3500 mm)
Execute moving Theta-D. (0.0000 deg)
Execute measuring Theta-S.
(START=-1,STOP=2.96,STEP=0.02,SPEED=4)
  
```


**Measurement Parameters:**

Zero Position	: -1.5000 mm	PHA Window	: 100 div
Theta-D	: 0.0000 deg	Monochro 2Theta	: --- deg
Theta-S	: 0.0000 deg	Monochro Theta	: --- deg
Theta-X	: --- mm	Monochro ThetaK	: --- mm
HV	: 695 V	tIntegral	: 0.7220 micro s
PHA Baseline	: 50 div	tDifferential	: 0.1833 micro s




Now remove the center slit and put the long range slit in filter box.


- Note the correct orientation (as shown in the figure on the right) for the long range slit.
- Confirm the pop-up dialog to start new measurement.
- The htm program will perform Theata D correction
- The htm program will perform 2 Theata correction



## Theta D alignment (correction)




## 2 Theta alignment (correction)







Now install the alignment jig (actually it is labelled as Setting jig). A metal plate with standard compressed silicon (Si) powder in the standard sample holder. Confirm OK to proceed with next measurement of Theta which should produce a symmetric triangle shown below.



- Now,
- Keep the setting jig
  - Remove the Al foil, if present (Absorber)
  - Remove the Ni foil, if present (Beta filter)
  - Remove the high limiting jig (03. HL)
  - Confirm OK to proceed to next measurement
  - That will be final alignment

**Intensity (cps)**

25000  
20000  
15000  
10000  
5000  
0

46.500 47.000 47.500 48.000

**2Theta/Theta (deg.)**

Grid Print Save Result Process Help Execute Counting Loss

2Theta Alignment Voltage 1: 40 kV Current 1: 40 mA HV/PHA Mono Sample d: 3.13532 A XG ending  
 Theta Alignment Voltage: HV/PHA Adjustment Monochro Sample d: 3.3756 A Present  
 Monochromator Adjustment Profile Sample d: 1.91999 A Aging  
 Profile Measurement HV: PHA Window: 80 div Xray off  
 Power off

Goniometer : Ultima3 theta-theta goniometer  
Attachment : Standard sample holder  
Slit : Auto Variable slit(Parallel/Ultima3)  
Counter : Scintillation counter  
Monochromator : Not installed  
X-ray tube : Cu (Line focus)

Zero Position:-1.1750 mm

Theta-X :  
HV :  
PHA Baseline :

-1.1900 mm PHA Window : 80 div  
-0.2523 deg. Monochro 2Theta: ... deg.  
1.1642 deg. Monochro Theta: ... deg.  
... mm Monochro ThetaX: ... mm  
695 V tIntegral : 0.7220 micro s  
60 div tDifferential : 0.1833 micro s

Execute setting XG power. (40 kV, 40 mA) My Computer

**Right Console**

Standby now !!

Pause Break

## COMPLETE ALIGNMENT, SAVE & PRINT RESULTS

At this point the alignment is complete. However, it is very important to “Save” the results. Now the “Save” and “Print” buttons will become active. First, click on “Save” button. This saves the current alignment values.


Now scroll down the .htm file’s left hand section and note the following values:

- Ts
- Zero Position
- Theta D Correction
- Theta S Correction

Now click on the “Print” button and select the “Registered Alignment Results”. The printer on top of the cabinet will print a page similar to the one shown on next page. Write down the four values we got above on the sheet of paper and save it in the drawer.



## AUTOMATIC ALIGNMENT (Right)




## SWAP THE EQUIPMENT BACK TO ORIGINAL CONDITION




Now having completed the alignment using the Scintillation Detector, we need to switch back to the original set-up with Mercury CCD.

- First, we need to bring the XG operation back to 20 kV and 2 mA. This can be done using the XG Operation Program
- Now we need shut-down the X-Ray Generator by doing the following:
  - Press the grey “OFF” button below X-RAYS in the front panel of the generator
  - Next press the red “OFF” button below POWER in the front panel of the generator
  - Switch “OFF” Servo AMP Unit after opening the bottom left panel
- Shut-down the Control PC
- Remove the Scintillation Detector by doing the following:
  - Un-plug the Scintillation Detector cable
  - Un-plug the RJ1 cable
  - While holding the Scintillation Detector, undo the screws that hold it and gently remove it
- While holding the Filter Holder, undo the screws that hold it and gently remove it and remove all the components inside it. Store the components inside the cabinet above the Control PC.
- Now, remove the two (2) black pieces that were holding the Filter Holder and the Scintillation Detector and place them along with their screws inside the cabinet above the Control PC.
- Now, remove the Standard Sample Holder by doing the following:
  - Remove the three thumb screws that hold the Sample Holder
  - Gently wiggle the Holder and pull toward you until comes out
  - Put the Holder in its holding-base and store it inside the cabinet above the Control PC
- Next, attach the black piece that holds the CCD on 2-Theta arm and secure all six screws
- Next, while carefully holding the CCD place it in the black stripe’s dove-tail and move the CCD until the mark on the top CCD matches the original distance before removing the CCD (in our case it was 105 mm)
- Now using the four Philips screws, while still holding the CCD with your one hand, fully secure the CCD the black stripe.
- Next, put the original MicroArea holder back in the middle of the goniometer (from where we just removed the Standard Sample Holder) by gently wiggling it inside the center hole
- Using the two hex screws secure the MicroArea holder ensuring the video camera is not pushed back or forward
- Swap the H10 slit with the original H2 in the CBO box
- Remove the CCD cover
- Close the X-Ray Cabinet door
- Switch the SERVO AMP unit power
- Now we need start the X-Ray Generator by doing the following:
  - Press the green “ON” button below POWER in the front panel of the generator
  - Wait for the READY green LED to light below the X-RAYS panel
  - Next, press the white “ON” button below X-RAYS in the front panel of the generator
  - Make sure, the EXT button on Shutter 2 is ON
- Now restart the Control PC
  - Log-in to Control PC
- Start the Frame Grabber PC
  - Log-in to Frame Grabber PC



- Wait for the DOS window to come on
- Switch back to Control PC
- In Control PC do the following:
  - Right Click RGHEExpand
  - Under RINT2200, select, "Show Window"
    - This will open a RINT2200 control window
  - Under Right, select, "Show Window"
    - This will open a Right control window
  - Click open Rigaku Control Panel Directory
  - In System Details do the following:
    - Switch from XG Server Mode: Standard to Std (Shutter relay)
    - Click OK to confirm it
    - The Right Control window will update
  - Now click open the RINT2200 button in the Rigaku Control Panel Directory
  - Now do the following switch backs:

<p>Now click on the RINT2200 Right System (but not the <del>RINT2200 Right System Add/Delete button</del>). This should open a new dialog box. In the following three tabs RINT2200 Right System Property need to be changed.</p> <ul style="list-style-type: none"> <li>● System Construction</li> <li>● X-Ray Beam Type</li> <li>● Geometry System</li> </ul>	
System Construction	
CHANGE FROM	CHANGE TO
<i>Goniometer: Ultima 3 theta-theta goniometer</i>	<i>Goniometer: Ultima 3 theta-theta goniometer+Mercury CCD</i>
<i>Attachement: Standard Sample Holder</i>	<i>Attachment: Micro area</i>
<i>Filter: K-beta filter</i>	<i>Filter: Not installed</i>
<i>Detector: Scintillation Counter</i>	<i>Detector: Mercury CCD</i>
<i>Temp Contoller:Not installed</i>	<i>Temp Contoller:Not installed</i>
<i>Incident Monochromateor: use</i>	<i>Incident Monochromateor: use</i>
<i>Slit: Auto Variable Slit (Parallel/Ultima3)</i>	<i>Slit: Auto Variable (Mercury CCD/general)</i>
<i>Monochromator: Not installed</i>	<i>Monochromator: Not installed</i>
X-Ray Beam Type	
CHANGE FROM	CHANGE TO
No change. <i>Focus Selection: Line; Wave length selection: K-alpha</i>	
Geometry System	
CHANGE FROM	CHANGE TO
<i>Geometry Information: Parallel Beam</i>	<i>Geometry Information: Micro Area</i>
After the change the three tabs should like the following	



- Click OK
- Click Yes to confirm it

## SUMMARY

This procedure is what is being followed at Florida State University Chemical Sciences Laboratory Powder X-Ray Diffraction Facility's Rigaku Ultima III w/ Mercury CCD. It is likely that your procedure is slightly different from mine. If you have questions, please contact me at [tsomasundaram@fsu.edu](mailto:tsomasundaram@fsu.edu) or at 1-850-644-6448.

P.S. The procedure is still beta version 0.9.6 and needs minor tweaks especially toward the last stage of restarting the CCD. I will update this procedure in early 2015 after redoing the alignment couple of more times and taking precise notes.

© 2013-14 Thayumanasamy Somasundaram | Kasha Laboratory

Version 0.9.6 dated November 18, 2014

Institute of Molecular Biophysics | Florida State University, Tallahassee, FL 32306-4380

[tsomasundaram@fsu.edu](mailto:tsomasundaram@fsu.edu) | [www.biophysics.fsu.edu/soma](http://www.biophysics.fsu.edu/soma) | 1-850-644-6448